

Enhance Your Pool Access Today!

aQuaTrek™

- Ladders
- Steps
- Ramps
- Wheel Chairs

"It's an excellent product and we love it! You do a great job and we wouldn't get nearly the use out of the pool that we do if we didn't have it. If you ever need an endorsement, I'll give one."

— Otis Reeves
Prattville YMCA

- Practical Solution for Aquatic Access
- Immitated But NEVER Duplicated
- Affordable Pricing
- Personalized Customer Service
- Custom Built for Every Pool
- Handcrafted Workmanship
- All Deck Customizations Included

Contact Us Or Your Favorite Pool Products Dealer!

Find Out How

YOUR POOL ACCESS CAN BE IMPROVED!

"We love them! (an AQ-7000 AquaTrek 7-Treaded Step and an AQ-350-WC AquaTrek Wheel Chair) We are having a state-wide aquatics event for Special Olympics because of these! Thank you!"

— Paula Ann Sorensen of Iowa Lakes Community College

CALL TODAY AND GET YOUR POOL ACCESS UP-TO-DATE!

"We love your AquaTrek Pool Access System, very much. Everyone uses it, even people who thought they'd never use it because it was for handicapped use. It's just so safe and easy. It's probably the best thing we've ever done for our pool. If there is ever anything we can do or say for you, it would be our pleasure. Your service and product are absolutely excellent. We'd be happy to talk to people and provide product recommendations."

— Mike Pierce of Point East Condominium

REHAB SYSTEMS, LLC
Building Excellence One Step At A Time

aQuaTrek™
Pool Access Systems

Aquatic therapy - A wave of the future that will only increase the need for practical and dignified pool access. AquaTrek is the solution!

Can your facility effectively meet it's client's needs without the AquaTrek Pool Access System?



Website: www.rehabsystems.net Toll Free: 800-726-8620
ORDER YOUR ACCESS SYSTEM TODAY!

AQUATREK POOL ACCESS SYSTEMS

- Durable & Lightweight
- Slip Resistant Step Surface
- Portable & Removable
- Custom Built For Each Pool
- Full, One-Year Warranty
- Noncorrosive Construction
- Customizable Riser Heights
- Walk Forward To Enter & Exit
- Safe & Easy To Use

The AquaTrek Ladders, Steps And Ramps

Do you want to create safe and practical access to your pool? People of all ages and physical abilities can enter and exit the pool with safety, ease and confidence while using the AquaTrek Ladders, Steps or Ramps.


The AquaTrek products are all easily removable for automatic pool cleaners, retractable pool covers and solar blankets. Wheels built right into the units allow staff members to roll the units away in minutes.

Wide treads and handrails provide the needed stability. Custom built construction of the product includes a noncorrosive metal skeleton encased in a rigid PVC skin. All hardware is 304 stainless steel.

Though imitated, the AquaTrek Ladders, Steps and Ramps have never been duplicated for product quality.

“It (an AQ-5600 AquaTrek Ladder) permits people to get into our pool WITHOUT using the lift chair. This is a confidence builder to people who are otherwise losing mobility”

– Jamie Endresu
Mon Valley YMCA



Above: The AquaTrek Ladder

AquaTrek Ladder and Step models vary depending on pool depth and desired riser height (distance between steps).



Above: The AquaTrek Ramp

AquaTrek Ladder available in 4, 5 or 6 step models.

AquaTrek Step available in 6, 7 or 8 step models.

AquaTrek Ramp available in standard or long.



Above: The AquaTrek Step

All AquaTrek Ladders, Steps and Ramps are customized to your pool deck and gutter system.



The AquaTrek Wheel Chair

- 7 and 24 Inch Mag Wheels
- Totally Submersible
- Footrest & Swing Arm Standard
- Safe & Comfortable
- Stainless Steel Axles & Fasteners
- Stainless Steel Wheel Rims Standard

“The AquaTrek Wheel Chair is lightweight, comfortable for patrons and easy to maneuver. Thanks for a good experience.”

– Diane Novotny
Jacksonville Parks & Rec

The AquaTrek Wheel Chair is the ultimate way to give your pool users assistance around and into the pool. It is totally submersible without becoming damaged by the water, making it an ideal aquatic access product.

Website: www.rehabsystems.net Toll Free: 800-726-8620

ORDER YOUR ACCESS SYSTEM TODAY!