




Let's walk together through the access improvement of a walk-in system on step at a time. Please bookmark our website, www.rehabsystems.net, as a resource for further details. Our website includes detailed descriptions, frequently asked questions and testimonials.


AquaTrek Ladder

*The AQ-5000
 5-Tread Ladder*




The AquaTrek Ladder is the dignified, forward-walking entry and exit to the pool that makes pool access safe and easy. Each AquaTrek Ladder is comfortable to use with wide treads and strong, metal-reinforced siderails. Stability is guaranteed due to deck mounting features, while portability is assured by hitchpins that separate the unit from the deck mounting position. Unlike other walk-in pool entry systems, the AquaTrek Ladder offers small, customizable distances between steps to guarantee a therapeutic system. To be therapeutic, a walk-in system needs to have the smallest possible distance between steps. The AquaTrek Ladder construction combines a

rigid, metal skeleton with a tough, furniture grade PVC skin. Each step tread includes a two-tone, nonskid surface for safety. Built-in wheels on the AquaTrek Ladder allow ease of movement of the unit in and out of the pool and around the deck area. The AquaTrek Ladder is available in a standard unit rated for 600-lbs user weight capacity, or a heavy-duty unit to extend the standard rating. Each AquaTrek Ladder has a one-year warranty backed by a manufacturer established for over 15 years. The AquaTrek Ladder is available in four, five or six step models.



AquaTrek Step

*The AQ-6000
 6-Tread Step*




The AquaTrek Step is the next level for walk-in access to your pool. Designed as a successor unit to the AquaTrek Ladder, the AquaTrek Step is the solution for deep pool situations that demand six treads, or more, to be therapeutic. It is a small distance between steps that defines a walk-in system as therapeutic. The AquaTrek Step is available in six, seven or eight tread models.



AquaTrek Ramp

The AQ-9000 Ramp




The AquaTrek Ramp is the practical solution for ramped access to a swimming pool. The AquaTrek Ramp is available in either a standard or long length. By using the AquaTrek Ramp systems, clients receive a gentle, sloped entry where they can either walk or be wheeled into the pool. The two versions of the AquaTrek Ramp available are offered for the same value pricing. The AQ-9000 AquaTrek Ramp Standard is the ramp meant for pools less than 47-inches deep from pool floor to deck surface. The AQ-9000L AquaTrek Ramp Long is a longer length AquaTrek Ramp meant for pools where the depth from pool floor to deck surface is 48-inches or

greater. The AquaTrek Ramp is designed with cross-bars to directly connect the bottom of the two siderails, adding an extra layer of strength and durability to the unit. A 30-inch tread width on the AquaTrek Ramp assures stability, unlike wider systems where 36-inch distances between the siderails allow for more potential instability. Elevation kits are also available to help reduce the slope of an AquaTrek Ramp system in a pool. The AquaTrek Wheel Chair is the perfect companion item for use with the AquaTrek Ramp.



AquaTrek Wheel Chair

*The AQ-350-WC
 Wheel Chair*



The AquaTrek Wheel Chair is totally submersible and perfect for your aquatic wheel chair needs. It is designed specifically for use in either a zero depth pool or on a ramp system like the AquaTrek Ramp. The bushings of the AquaTrek Wheel Chair are custom designed using Nylatron to prevent corrosion. A retractable footrest and soft, cushioned seat on the AquaTrek Wheel Chair provide comfort for each client. Stainless steel handrims on the rear wheels allow client mobility. A swing-arm allows for ease of mobility on and off the AquaTrek Wheel Chair. The Furniture

grade PVC construction of the AquaTrek Wheel Chair ensures durability for the heavy-duty needs of the commercial aquatic environment. The options of a shoulder harness and seat belt can be purchased for the AquaTrek Wheel Chair upon request.

